DESCRIPTION DU DISPOSITIF

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DESCRIPTION OF THE DEVICE

The devices include a bag designed with one or two connection modes, according to the references. Each bag contains a certain volume of platelet additive solution SSP or SSP+ according to following composition:

- Composition SSP:** Na₂citrate 2H₂O: 2.94g, Naacetate 3H₂O: 4.08g, NaCl: 6.75g, Water for Injection 1000 ml, pH 7.2
- Composition SSP+:** Na₂citrate 2H₂O: 3.18g, Naacetate 3H₂O: 4.42g, NaH₂PO₄ 2H₂O: 1.05g, Na₂HPO₄ 3.05g, KCl: 0.37g, MgCl₂ 6H₂O: 0.30g, NaCl: 4.05g, Water for Injection 1000 ml, pH 7.2

The bag containing the solutions is non-PVC. The PVC of the tubings is DEHP plasticized. The biological products are not-PVC or circulating in the device could prevent container des phalates.

INDICATIONS

The platelet additive solutions SSP and SSP+ are designed to partially replace plasma in the preparation and storage of buffy-coat derived platelet concentrates or apheresis platelet units. The recommended ratios are as follows: up to 70% SSP / 30% plasma, up to 80% SSP+ / 20% plasma. These solutions enable platelets to be stored at 22 °C ± 2 °C, under gentle agitation, for up to 5 days (SSP) or 7 days (SSP+) following collection and according to local regulations.

SIDE EFFECTS

There are no known side effects with the use of SSP or SSP+.

WARNING

- Verify the integrity of the device and the overwrap before use.
- Use the device at a temperature between 18° and 24 °C
- Do not store the device at a temperature above 35 °C
- Do not use sharp objects.
- Do not use if the overwrap or the device show any visible signs of deterioration.
- Do not use if the solution is cloudy.
- Do not use for the transfusion of platelets.
- Do not directly infuse the solution.
- The device must be used with other approved devices including connectors complying with current standards: ISO 5942.
- Do not reuse the device, as there could be a risk of contamination. And additionally, the remaining quantity of solution might be insufficient.

INSTRUCTIONS FOR USE

- Manifold pooling (use of 4 to 6 Buffy-Coats)**
 - Make a sterile connection between each buffy-coat tubing and each end-line tubing of the manifold device.
 - Connect by sterile connection the SSP or SSP+ bag tubing to the remaining manifold set end-tubing.
 - Hang all the buffy-coat bags at the same level and the SSP or SSP+ bag on a higher level.
 - Allow the buffy-coats to drain into the pooling bag.
 - Close the clamp above the pooling bag and then allow the SSP or SSP+ solution to drain into each buffy-coat.
 - Agitate in order to recover a maximum of buffy-coat & open the clamp, allow the SSP/SSP+ solution to drain into the pooling bag.
 - Remove the manifold set and then follow standard procedures to obtain a pooled buffy-coat derived platelet concentrate.

Cascade pooling (4-6 Buffy Coats)

- Connect by sterile connection the SSP or SSP+ bag tubing to the bleed line of the first buffy-coat bag in the train.
- Make a sterile connection between the buffy-coat plasma line (tube equipped with a breakaway cannula) of the first buffy-coat bag and the bleed line of the second buffy-coat.
- Then make a sterile connection between the second buffy-coat and the third buffy-coat. Repeat the operation with each buffy-coat until a maximum of 6 buffy-coats are pooled.
- Operate the tubing of the plasma line of the last buffy-coat bag to a platelet recovery set.
- Hang the chain vertically by suspending the SSP or SSP+ bag with the chain of buffy-coats underneath. Allow the buffy-coats to drain into the pooling bag.
- Use the SSP or SSP+ solution to drain one bag after each other.
- Remove the set and then follow standard procedures to obtain a pooled buffy-coat derived platelet concentrate.

Apheresis

- Connected to the apheresis set, the SSP or SSP+ solution can be used as a substitution solution for a part of the plasma in platelet concentrates.
- Connect the Additive Solution with the luer connector to the apheresis set.
- Break the breakable connector in the Additive Solution line by moving the connector back and forth 90° until an audible "snaps" sound is heard in both directions (see drawing (A)).
- Band tubing 180° to force the arrow piece of the frangible upward into the soft tubing closer to the bag opening and away from the base (B). This prevents this piece from resealing. To ensure proper flow, check to see if a gap exists between the arrow and base of the frangible (C). If not repeat step B.
- Hang the Additive Solution bag on the IV hook with the platelet product bag supported (not hanging) below.
- Follow the valid apheresis protocol of your institution. Allow the appropriate amount of storage solution to flow into the platelet product bag and mix gently.

DISCARD WARNINGS

Discard according to the local regulation(s) on the potentially dangerous biological products.

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