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Anticoagulant solutions for automated apheresis

ENGLISH
- Sterile and non pyrogenic anticoagulant solutions for blood anticoagulation.
- MacoPharma anticoagulant solutions are strictly reserved for the automated apheresis procedures.

✓ **COMPOSITION OF THE BAG**
- Plastified PVC bag.

✓ **COMPOSITION OF ACD formula A**
Trisodium Citrate, 2H₂O: 22.0g, Citric Acid, 1H₂O: 8.0g, Glucose, 1 H₂O: 24.5g, PPI Water: 1000ml

✓ **COMPOSITION OF CITRATE 4%**
Citric Acid, 1H₂O: 2.8g, Trisodium Citrate, 2H₂O: 39.6g, PPI Water: 1000ml

✓ **PRECAUTIONS FOR USE**
- Do not use for direct infusion.
- Do not use if solution is cloudy, or if there is any visible sign of deterioration.
- Do not re-use the bag. After opening the bag, discard the non used product.
- Do not vent.
- Do not use if the overwrap and/or bag are damaged.
- Check expiry date before use.
- Septic risk, risk of misuse.
- The PVC complies with the European Pharmacopoeia. The plasticizer is DEHP. The biological products stored and/or circulating in these devices could contain phthalates.
- Reject according to local regulations' on the potentially dangerous biological products

✓ **STORAGE CONDITIONS**
- 2 years shelf-life in original sealed overwrap.
- Do not store above +35°C in original overwrap.
- Do not freeze.
- Protect from light during storage.

✓ **DOSAGE**
Anticoagulant solutions are added to whole blood during collection in a proportion defined by the person responsible for the procedure. The following rates are commonly advised by apheresis machines manufacturers : 1/9 for ACD-A solution, 1/15 for citrate 4%. This rate, however, can be changed according to the performance of the apheresis machines.

✓ **Instructions for use**
Installation of the apheresis kit on the separator: the connection of the MacoPharma anticoagulant solution bag to the apheresis system is carried out according to the separator's instructions, by 2 means :
- break and detach the twist-off situated at the end of the anticoagulant solution, and connect the bag using a spike.
- connect the male luer, located on the anticoagulant solution.
If you use the luer connection:

1. Connect the Anticoagulant Solution with the luer connector to the apheresis set.
2. Break the breakable connector in the Anticoagulant Solution line by moving the connector back and forth 90° until an audible "snap" sound is heard in both directions. This separates the top piece of the frangible (arrow-shaped design) from the male luer base.
3. Bend tubing 180° to force the arrow piece of the frangible upward into the soft tubing closer to the bag opening and away from the base. This prevents this piece from reseating. To ensure proper flow, check to see if a gap exists between the arrow and base of the frangible. If not repeat step 3.
4. Hang the Anticoagulant Solution bag on the IV hook.
5. Follow the valid apheresis protocol of your institution.

✓ **CONTRA-INDICATIONS**
Contra-indications are those of blood collection. For cellular apheresis, check haemostasis. Use only for automated plasmapheresis or cytopheresis procedures. Do not perform apheresis in case of: serious haemopathy with haemostatic trouble, cerebral vascular accidents, pericardial or endocardial cardiopathy, post surgical period after intervention on central nervous system and spinal cord.

✓ **UNDESIRABLE EFFECTS**
Paresthesia : stinging sensations on lips and extremities ; Nausea ; Dysgueusia : weird taste sensation ; Cold sensation.

✓ **PREGNANCY AND LACTATION**
Do not carry out apheresis among pregnant women, nor during lactation (except plasmatic exchange if absolutely (indicated) necessary).

✓ **OVERDOSE**
Tetany possibility ; vomiting ; hypotension ; tachycardia. In case of overdose, control clinic and electrographic signs. If necessary inject calcium until (correction of the troubles) symptoms disappear.

✓ **SPECIAL WARNINGS**
Calcium injection is often necessary during cellular apheresis to limit the drop of the ionised calcium rate.

✓ **EFFECTS ON ABILITY TO DRIVE**
Rest time advised after an apheresis collection.

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