
MACOPHARMA
Platelets Solution

THERAFLEX UV-Platelets

Light the way to blood safety

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
THERAFLEX UV-Platelets:

Leistung

Viren ^(6 – 16)	Log-Reduktion
Murine Leukemia virus (MuLV)	< 1
Human immunodeficiency virus (HIV)	1.0
Bovine viral diarrhea (BVDV)	≥ 2.7
Suid herpesvirus 1 (SHV-1)	2.8-3.9
Feline calici virus (FCV)	3.0 ± 0.2
Crimean-Congo hemorrhagic fever virus (CCHFV)	≥ 2.2
West Nile Virus (WNV)	3.6 - 4.4
Severe Acute Respiratory Syndrom virus (SARS-CoV)	≥ 3.4
Middle East respiratory syndrome coronavirus (MERS-CoV)	≥ 3.7
Hepatitis E virus (HEV)	≥ 3.5
Encephalomyocarditis virus (EMCV)	3.0-3.7
Hepatitis A virus (HAV)	≥ 4.2
Dengue virus 2-4 (DENV)	4.4-5.4
Ebola virus (EBOV)	≥ 4.5
Usutu virus (USUV)	5.2 ± 0.7

Viren ^(6 – 16)	Log-Reduktion
Nipah virus (NiV)	≥ 4.3
Yellow fever virus (YFV)	≥ 4.8
Porcine Parvo virus (PPV)	≥ 6.1
Hepatitis C virus (HCV)	≥ 5.0
ZIKA virus (ZIKV)	5.7
Severe Acute Respiratory Syndrom virus-2 (SARS-CoV-2)	≥ 5.2
Ross River virus (RRV)	> 5.1
Sindbis virus (SIV)	≥ 5.2
Influenza virus (H3N2)	≥ 5.3
Dengue virus 1 (DENV)	≥ 5.3
Vesicular stomatitis virus (VSV)	≥ 6.3
Chikungunya virus (CHIKV)	6.3
Japanese encephalitis virus (JEV)	7.0

Parasiten ¹⁸	Log-Reduktion
<i>L. infantum</i>	> 1.4
<i>T. cruzi</i> (Chagas Disease)	> 2.8- > 4.2
<i>B. divergens</i>	4 -5
<i>P. falciparum</i>	> 5

THERAFLEX UV-Platelets:

Leistung

Bakterien ⁽¹⁷⁾	Log-Reduktion	Bakterien ⁽¹⁷⁾	Log-Reduktion
<i>Enterobacter cloacae</i> *	6.3 ± 0.6	<i>Streptococcus bovis</i> *	7.0 ± 0.3
<i>Escherichia coli</i> *	7.3 ± 0.4	<i>Streptococcus dysgalactiae</i> *	3.5 ± 0.3
<i>Klebsiella pneumoniae</i> *	5.9 ± 0.7	<i>Streptococcus pyogenes</i> *	4.4 ± 0.8
<i>Morganella morganii</i> *	7.5 ± 0.3	<i>Listeria monocytogenes</i>	4.8 ± 0.3
<i>Proteus mirabilis</i> *	7.0 ± 0.3	<i>Acinetobacter baumannii</i>	5.1 ± 0.3
<i>Pseudomonas fluorescens</i> *	7.0 ± 0.2	<i>Streptococcus agalactiae</i>	≥ 6.0
<i>Serratia marcescens</i> *	5.8 ± 0.2	<i>Streptococcus pneumoniae</i>	4.8 ± 1.0
<i>Staphylococcus aureus</i> *	6.6 ± 0.4	<i>Bacillus cereus</i> *	3.1 ± 0.1
<i>Staphylococcus epidermidis</i> *	4.6 ± 0.5	<i>Bacillus thuringiensis</i> *	3.6 ± 0.4
<i>Clostridium perfringens</i> *	≥ 4.7	<i>Propionibacterium acnes</i>	4.5 ± 1.1
<i>Pseudomonas. Aeruginosa</i>	≥ 4.9		

Leukozyten ⁽¹⁹⁾	Log-Reduktion
In vitro (Proliferation/ Limiting dilution assay)	> 5.0
In vivo (Mouse Model)	GvHD verhindert

Log-Reduktion

*WHO International Repository Platelet Transfusion Relevant Reference Strains

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
THERAFLEX
Platelets
LIGHT THE WAY TO BLOOD SAFETY

“*Solutions to make
the best out of every
drop of blood.*”

Literaturverzeichnis :

1. Mohr H, Gravemann U, Bayer A, Müller TH: Sterilization of platelet concentrates at production scale by irradiation with short-wave ultraviolet light. *Transfusion* 2009;49:1956–1963.
2. Mohr H, Gravemann U, Müller TH: Inactivation of pathogens in single units of therapeutic fresh plasma by irradiation with ultraviolet light. *Transfusion* 2009;49:2144–2151.
3. Mohr H, Steil L, Gravemann U, Thiele T, Hammer E, Greinacher A, Müller TH, Volker U: A novel approach to pathogen reduction in platelet concentrates using short-wave ultraviolet light. *Transfusion* 2009;49:2612–2624.
4. Jerard Seghatchian, Frank Tolkdorf. Characteristics of the THERAFLEX UV-Platelets pathogen inactivation system – An update. *Transfusion and Apheresis Science* 46 (2012) 221–229.
5. Data on file at Macopharma.
6. Mohr, H., et al. (2009). «A novel approach to pathogen reduction in platelet concentrates using short-wave ultraviolet light.» *Transfusion* 49(12): 2612-2624.
7. Steinmann, E., et al. (2013). «Two pathogen reduction technologies--methylene blue plus light and shortwave ultraviolet light--effectively inactivate hepatitis C virus in blood products.» *Transfusion* 53(5): 1010-1018.
8. Faddy, H. M., et al. (2016). «Inactivation of dengue, chikungunya, and Ross River viruses in platelet concentrates after treatment with ultraviolet C light.» *Transfusion* 56(6 Pt 2): 1548-1555.
9. Fryk, J. J., et al. (2017). «Reduction of Zika virus infectivity in platelet concentrates after treatment with ultraviolet C light and in plasma after treatment with methylene blue and visible light.» *Transfusion* 57(11): 2677-2682.
10. Eickmann, M., et al. (2018). «Inactivation of Ebola virus and Middle East respiratory syndrome coronavirus in platelet concentrates and plasma by ultraviolet C light and methylene blue plus visible light, respect.»
11. Faddy, H. M., et al. (2019). «Inactivation of yellow fever virus in plasma after treatment with methylene blue and visible light and in platelet concentrates following treatment with ultraviolet C light.» *Transfusion* 59(7): 2223-2227.
12. U., et al. (2018). «Ultraviolet C light efficiently inactivates nonenveloped hepatitis A virus and feline calicivirus in platelet concentrates.» *Transfusion* 58(11): 2669-2674.
13. Eickmann, M., et al. (2020). «Inactivation of three emerging viruses - severe acute respiratory syndrome coronavirus, Crimean-Congo haemorrhagic fever virus and Nipah virus - in platelet concentrates by ultraviolet C light and in plasma by methylene blue plus visible light.» *Vox Sang.*
14. Praditya, D., et al. (2020). «Hepatitis E virus is effectively inactivated in platelet concentrates by ultraviolet C light.» *Vox Sang* 115(7): 555-561.
15. Schulze, T. J., et al. (2022). «THERAFLEX ultraviolet C (UVC)-based pathogen reduction technology for bacterial inactivation in blood components: advantages and limitations.» *Annals of Blood*; Vol 7 (September 30, 2022).
16. Rustanti, L., et al. (2020). «Inactivation of Japanese encephalitis virus in plasma by methylene blue combined with visible light and in platelet concentrates by ultraviolet.»
17. Gravemann U, Handke W, Müller TH, Seltsam A. Bacterial inactivation of platelet concentrates with the THERAFLEX UV-Platelets pathogen inactivation system. *Transfusion*. 2019 Apr;59(4):1324-1332
18. Castro, E., et al. (2014). «The efficacy of the ultraviolet C pathogen inactivation system in the reduction of *Babesia divergens* in pooled buffy coat platelets.» *Transfusion* 54(9): 2207-2216.
19. Pohler, P., et al. (2015). «Pathogen reduction by ultraviolet C light effectively inactivates human white blood cells in platelet products.» *Transfusion* 55(2): 337-347.

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