




## COLLECTION LINE



# COLLECTION LINE OPTIMAL COLLECTION PROCEDURE

## Permanent Clamps (a)

- Avoid bacterial contamination
- Replaces need for Grommet
- Easy to use (one hand)

## Vacuvam (b)

- Sampling with vacuum tubes
- Closed system
- Hygienic & safe
- Integrated device
- Easy to use

## Bactivam (c)

- Diversion of first mL of blood
- Guaranteed Sample availability
- Reduction of bacterial contamination
- Optimized sample collection
- Full traceability


## 16-Gauge Needle (d)

- 100% quality controlled for sharpness
- Double silicon coated, angled bevel for smooth insertion.
- Ergonomic hub
- Lockable in Secuvam for full safety

## Secuvam (e)

- Prevention of needle stick injury
- Round shape for ease of use
- Integrated device

**FULLY  
INTEGRATED  
SAFE AND EASY  
TO USE**





Distributed by MACOPHARMA USA  
3675 Crestwood Parkway, Suite 260  
Duluth, GA 30096

Manufactured by MACOPRODUCTIONS SAS  
200, chaussée Fernand Forest  
59200 Tourcoing - France

**TRANSFUSION**  
**macopharma**  
DESIGNED FOR LIFE