

Macopharma Press Release

Monday April 27, 2020

Macopharma announces implementation of Correct Connect System

Macopharma is pleased to inform that « Correct Connect » anticoagulant solution ACD-A used on apheresis devices recently acquired the CE marking by its notified body so called G-MED, in accordance to ISO 18250-8 norm.

ISO 18250-8 norm was published in August 2018, coming out from an « apheresis industry working group » aiming to secure apheresis collections.

The idea of developing standardized connection solutions for the apheresis industry was initiated in 2014 and the working group was led by Terumo BCT. Participating members include Cerus, the European Blood Alliance (EBA), Fresenius Kabi, Grifols, Haemonetics, Macopharma, Qosina, Scinomed, Technoflex, Terumo BCT and Therakos.

The Correct Connect connector is a female connector that was designed in order to reduce the risk of misconnections and secure donor's safety.



The new ACD-A bag systems are easily identified by a specific labelling including the « Correct Connect » pictogram. The use of Correct Connect ACD-A bag system requires an adaptable apheresis disposable kit to be connected to, meeting as well the requirements of the ISO 18250-8 norm.



We, as Macopharma, stay at your disposal for any further question.